Giardia Single Test 10 x #(2)800009
Rev.: KI-40-00510-01

Giardia Single Test

Rapid Test

- (PT) Teste de Giardia
- (CS) Giardia Single Test
- (PL) Test Giardia
- (RO) Giardia Test unic
- (SK) Giardia Single Test
- (DA) Giardia Single Test
- (NO) Giardia enkelttest
- (SV) Giardia test, enpack
- (ZH) 贾第鞭毛虫单项测试

Giardia Single Test

Rapid Test

- (ES) Test único de Giardia
- (FR) Giardia Single Test
- (DE) Giardia Einzeltest
- (IT) Giardia Single Test
- (NL) Giardia enkele test

For use in dogs and cats

Reorder # (2)800009

LOT 40005TTMMJJ

2022-03


30°C
2°C

8 720171 399919

Distributed by:
Covetrus BV
Beversestraat 23
5431 SL Cuijk (NL)
cbproducts@covetrus.com
Distrivet, S.L.
Registro entidad
importadora HCMR-0097

Giardia Single Test 2 x #(2)800019
Rev.: KI-40-00502-01


Giardia Single Test

Rapid Test

- (ES) Test único de Giardia
- (FR) Giardia Single Test
- (DE) Giardia Einzeltest
- (IT) Giardia Single Test
- (NL) Giardia enkele test
- (PT) Teste de Giardia
- (CS) Giardia Single Test
- (PL) Test Giardia
- (RO) Giardia Test unic
- (SK) Giardia Single Test
- (DA) Giardia Single Test
- (NO) Giardia enkelttest
- (SV) Giardia test, enpack
- (ZH) 贾第鞭毛虫单项测试

For use in dogs and cats

Reorder # (2)800019

LOT 40005TTMMJJ

2022-03

8 720171 399919

Distributed by:
Covetrus BV
Beversestraat 23
5431 SL Cuijk (NL)
cbproducts@covetrus.com
Distrivet, S.L.
Registro entidad
importadora HCMR-0097

Rev.: PO-40-005-01

Giardia Single Test single teststrip

LOT: 40005TTMMJJ Exp.: 2022-03

Distributed by:
Covetrus BV
Beversestraat 23
5431 SL Cuijk (NL)
cbproducts@covetrus.com

Giardia Single Test

Instructions for use

- ES** Test único de Giardia · Instrucciones de uso
- FR** Giardia Single Test · Mode d'emploi
- DE** Giardia Einzeltest · Gebrauchsanweisung
- IT** Giardia Single Test · Istruzioni per l'uso
- NL** Giardia enkele test · Gebruiksaanwijzing
- PT** Teste de Giardia · Instruções de uso
- CS** Giardia Single Test · Návod na použití

EN Giardia Single Test Instructions for use

The Giardia test is used to test rapidly and reliably for Giardia antigen in faecal samples.
Handling the Sample
 Fresh faecal requires no pre-treatment and can be examined directly. The samples can be cooled between 2-8 °C for up to two days and stored at -20 °C for several weeks. Before the test is carried out, the samples must be brought to room temperature. Frozen samples should be thoroughly mixed after thawing.

Test procedure (Figure 1-4)

- Figure 1: Take a faecal sample using the cotton swab. Make sure that the upper end of the cotton swab is covered.
- Figure 2: Unscrew the test tube and place the faecal cotton swab inside. Use gentle pressure to apply an initial three drops of sample material to the sample well of the test cassette. If the liquid does not run up into the sample well, add one more drop of the sample liquid into the sample well. If the liquid remains in the sample well and does not draw up, it helps to press slightly on the sample field with the back of the cotton swab.
- Figure 4: Break of the top of the test tube by applying strong pressure. Open the test cassette packaging and place the cassette on a horizontal surface. Use gentle pressure to apply an initial three drops of sample material to the sample well of the test cassette. If the liquid does not run up into the sample well, add one more drop of the sample liquid into the sample well. If the liquid remains in the sample well and does not draw up, it helps to press slightly on the sample field with the back of the cotton swab.

Analysis (Figure 5 and 6)
 After 10 minutes one or two blue lines will appear in the reaction field, and the test can be analysed.

Positive test result (Figure 5):
 The test line and the control line become visible. Giardia antigens are detected. NB: Even a faint test line should not be interpreted as a positive test result. An excessive amount of sample feces can cause a brown line in the T-region (T), which must not be seen as a positive result. The test is invalid and should be repeated.

Negative test result (Figure 6):
 Only the control line becomes visible.

Invalid:
 If the control line is shown, the test is invalid. NB: All test results that are read later than 10 minutes are invalid.

Important:
 Asyst discharge does not occur with every passing of faecal, a negative test result does not necessarily preclude a repeat test. The test should be carried out once more with a pooled faecal sample after approximately one week.

Disposal
 All materials used must be disposed of carefully, because the material could be infectious. Test materials and used cassettes must be put in plastic bags and then disposed of.

Please note
 For veterinary use only. • Only for professional use. • Use the cassette as soon as possible after purchase. • Please refer to the instructions for required sample quantity. • Do not apply any sample material to the reaction field. • Do not touch the reaction field. • Use only the fluid supplied. • To avoid cross-contamination, a separate collection vessel should be used for each sample. • Faecal samples may be infectious, so suitable precautions should be taken for handling and disposal: avoid skin contact, wear laboratory gloves and clothing. • Do not use after the expiry date. • Do not use the test cassette if the packaging is damaged. • Observe analysis times indicated. • Store and transport the test only at the temperatures indicated.

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