CPV + FCoV + Giardia Combi Test

Rev.: KI-40-10210-01

10 x #(2)800011

CPV + FCoV + Giardia Combi Test

Rapid Test

- (PT) Teste Combinando CPV + FCoV + Giardia
- (CS) Kombinovaný test CPV + FCoV + Giardia
- (PL) CPV + FCoV + Giardia Combi Test
- (RO) CPV + FCoV + Giardia Combi Test
- (SK) Kombinovaný test CPV + FCoV + Giardia
- (DA) CPV + FCoV + Giardia-kombitest
- (NO) CPV + FCoV + Giardia Combi Test
- (SV) Kombinationstest för CPV + FCoV + Giardia
- (ZH) CPV + FCoV + Giardia 联合测试

CPV + FCoV + Giardia Combi Test

Rapid Test

- (ES) Test combinado de PVC + CoVF + Giardia
- (FR) Test combiné CPV + FCoV + Giardia
- (DE) CPV + FCoV + Giardia Combi Test
- (IT) Test combinato CPV + FCoV + Giardia
- (NL) CPV + FCoV + Giardia -combitest


For use in dogs and cats

Reorder # (2)800011

LOT 40102TTMMJJ

2022-03

30°C
2°C


Distributed by:
Covetrus BV
Beversstraat 23
5431 SL Cuijk (NL)
cbproducts@covetrus.com
Distrivet, S.L.
Registro entidad
importadora HCMI-0097

CPV + FCoV + Giardia Combi Test


Rev.: KI-40-10202-01

2 x #(2)800021

LOT
40102TTMMJJ

2022-03

CPV + FCoV + Giardia Combi Test

Rapid Test

- (ES) Test combinado de PVC + CoVF + Giardia
- (FR) Test combiné CPV + FCoV + Giardia
- (DE) CPV + FCoV + Giardia Combi Test
- (IT) Test combinato CPV + FCoV + Giardia
- (NL) CPV + FCoV + Giardia-combitest
- (PT) Teste Combinando CPV + FCoV + Giardia
- (CS) Kombinovaný test CPV + FCoV + Giardia
- (PL) CPV + FCoV + Giardia Combi Test
- (RO) CPV + FCoV + Giardia Combi Test
- (SK) Kombinovaný test CPV + FCoV + Giardia
- (DA) CPV + FCoV + Giardia-kombitest
- (NO) CPV + FCoV + Giardia Combi Test
- (SV) Kombinationstest för CPV + FCoV + Giardia
- (ZH) CPV + FCoV + Giardia 联合测试

For use in dogs and cats Reorder # (2)800021

Rev.: PO-40-102-01

covetrus

CPV + FCoV + Giardia Combi Test triple teststrip

LOT: 40101TTMMJJ Exp.: 2022-03

Distributed by:
Covetrus BV
Beversstraat 23
5431 SL Cuijk (NL)
cbproducts@covetrus.com

CPV + FCoV + Giardia Combi Test**Instructions for use****(ES) Test combinado de PVC + CoVF + Giardia • Instrucciones de uso****(FR) Test combiné CPV + FCoV + Giardia • Mode d'emploi****(DE) CPV + FCoV + Giardia Combi Test • Gebrauchsanweisung****(IT) Test combinato CPV + FCoV + Giardia • Istruzioni per l'uso****(NL) CPV + FCoV + Giardia-combitest • Gebruiksaanwijzing****(PT) Teste Combinando CPV + FCoV + Giardia • Instruções para uso****(CS) Kombinovaný test CPV + FCoV + Giardia • Návod na použití****(EN) CPV + FCoV + Giardia combined test****Instructions for use**

The Covetrus CPV+FCoV+Giardia combined test is a fast and reliable detection of parvovirus (CPV) and FCoV/coronavirus (CCV and FCoV) and Giardia antigens in faecal samples.

Handling the Sample

Fresh faecal requires no pre-treatment and can be examined directly. The samples can be stored between 2-8°C for up to two days and at -20°C for several weeks. Before the test is carried out, must be brought to room temperature. Frozen samples should be thoroughly mixed after thawing.

Test procedure (Figure 1-4)

Figure 1: Take a faecal sample using the cotton swab. Make sure that the upper end of the cotton swab is covered.

Figure 2: Unscrew the test tube and place the faecal cotton swab with the faeces inside.

Figure 3: Move the cotton swab back and forth in the liquid so that the faeces can dissolve completely. After about 10 seconds, add one more drop of the sample liquid into the sample well.

NB: Please do not allow the sample to sit after it has been dissolved. Apply it immediately to the sample well.

Figure 4: Break the spike of the test tube by applying pressure sharply. Open the test cassette package and place the cassette on a horizontal surface. Use the dropper to add one drop of the sample liquid into the sample well. Add one sample well of the test cassette. If the liquid does not start to run up after some seconds, add one more drop of the sample liquid into the sample well. If the liquid remains in the sample well and does not draw up, it helps to press lightly on the sample well with the back of the cotton swab.

Analysis (Figure 5)

After 10 minutes, one or two red lines appear in the reaction field, and the test can be analysed.

Positive test result (Figure 5 - +):

The test line and the control line become visible, and parvovirus is detected. The figure shows a positive parvovirus test result.

NB: Even a faint test line must be interpreted as a positive test result.

An invalid result (Figure 5 - -):

Only the control line becomes visible. The figure shows a negative coronavirus test result.

An invalid result (Figure 5 - 0):

If no control line is shown, the test is invalid.

NB: All results that are read later than after 10 minutes are invalid.

Important:

As cyst discharge does not occur with every passing of fecal, a negative test result does not necessarily preclude a persistent Giardia duodenalis infection. The test should be carried out once more with a pooled faecal sample after approximately one week.

Disposition

All materials used must be disposed of carefully, because the material could be infectious. Test materials and used cassettes must be put in plastic bags and then disposed of.

Please note

For veterinary use only. • Only for professional use. • Use the cassette as soon as possible after opening. • Please follow the procedure for required sample quantity. • Do not apply any sample material to the reaction field. • Do not touch the reaction field. • Use only the fluid from the sample well. • To avoid cross-contamination, use a new cassette and a new dropper for each sample. • Stool samples may be infected, so suitable precautions should be taken for handling and disposal avoid skin contact and wear laboratory gloves and clothing. • Do not use after the expiry date. • Do not use the test cassette if the packaging is damaged. • Observe analysis times indicated. • Store and transport the test only at the temperatures indicated.

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